

BED

SOLUTIONS

BED WITH BACKREST ON RATCHET AND PLAIN BED



Mi-7002 CX

Bed with Backrest on Ratchet

Introduction:

Medimek is a trusted brand in the medical fraternity over a decade. Our singular and simple philosophy is to cater to the needs of all our customers with a comprehensive range of excellent quality, user friendly hospital furniture and allied products, offered in a broad spectrum of exciting shades, aesthetic patterns and creative designs and backed by prompt after sales services.

Beds are universal and essential for ensuring successful medical treatment to the patient in any hospital and medical institution.

No hospital infrastructure is complete without the **Plain Beds**; basic, conventional and economical beds that are strong, sturdy but have mechanically simple operations. Largely opted for in the rural hospitals in the remote and inaccessible areas, they are also offered in customized sizes.



Mi-7003

Bed with Backrest on Ratchet
with SS head and foot boards with SS tubes



Mi-7012 CX

Plain Bed with SS head and foot boards
with Metal panels and Collapsible railings



Mi-7013

Plain Bed

with SS head and
foot boards with
SS tubes



Mi-7041

Pediatric Bed

with MS bows
and Full length
drop side
railings



Backrest

Options Available :

Manually Operated Beds with Backrest on Ratchet:

- Backrest adjustment on ratchet mechanism with five inclinations
- Two section uniformly perforated mild steel CRCA sheet top

Plain Bed:

- Plain uniformly perforated mild steel CRCA sheet top

Common Features:

- Top frame made from mild steel tubes
- "H" type tubular legs fitted with rubber shoes
- Four locations on bed frame for IV rod
- Pretreated and powder coated finish
- Supplied in SKD condition

Bed dimensions :

Option 1:

- Bed frame approx. dimension : 1980 mm L x 915 mm W x 600 mm H **[78"L x 36"W]**
- Overall approx. dimension (buffer-to-buffer) : 2030 mm L x 960 mm W

Option 2:

- Bed frame approx. dimension : 1980 mm L x 840 mm W x 600 mm H **[78"L x 33"W]**
- Overall approx. dimension (buffer-to-buffer) : 2030 mm L x 890 mm W

Option 3:

- Bed frame approx. dimension : 1905 mm L x 840 mm W x 600 mm H **[75"L x 33"W]**
- Overall approx. dimension (buffer-to-buffer) : 1960 mm L x 890 mm W

Head and Foot boards:

- Stainless steel head and foot boards with Metal panels / Stainless steel head and foot boards with SS tubes

Pediatric Bed with MS bows and full length drop side railings:

- Overall approx. dimension: 1450 mm L x 880 mm W
- Bed frame approx. dimension: 1370 mm L x 750 mm W X 590 mm H
- Plain uniformly perforated mild steel CRCA sheet top
- Both head bow and leg bow of equal height of approx. 1090 mm H
- Full length drop type side railing (pair)
- Two locations for IV rod
- Legs fitted with rubber shoes
- Pretreated and powder coated finish

Optional Accessories:

Mi-1007

Collapsible railings with SS tubes and MS brackets (pair)

Mi-1008

Swing away type SS tubes railings (pair)

Mi-1008 A

Swing away type MS tubes railings (pair)

Mi-1022

Two section 100 mm thick PU foam mattress covered with rexine

Mi-1023

Plain 100 mm thick PU foam mattress covered with rexine

Mi-1026

Pillow covered with rexine

Mi-1067

Telescopic IV Pole with two hooks - SS

Mi-1068

IV rod with two hooks - SS

Mi-1059

Four heavy duty 125 mm dia castors, two with brake

Mi-6060

Crib with Attachment

Mi-1082

Chart Holder

- Accessories, colors and designs shown in the brochure illustrations may not be the part of standard product configuration. • We reserve the right to make changes to designs, specifications and illustrations without prior notice. • Reproduction of this brochure in any form is liable for legal action.