

Mi-9005 A

**Electric ICU Bed - Regular Model
(with ACP box)**



Mi-9005 B

**Electric ICU Bed - Regular Model
(without ACP box)**



TECHNICAL SPECIFICATION

Product Name and Description	Product Code	
	Mi-9005 A	Mi-9005 B
ELECTRIC ICU BED – REGULAR MODEL		
Fully Motorized height adjustable ICU Bed	✓	✓
Overall approximate dimension (buffer-to-buffer): 2200 mm L x 1000 mm W	✓	✓
Bed frame approximate dimension: 2105 mm L x 930 mm W.	✓	✓
Mattress approximate dimension: 1980 mm L x 885 mm W.	✓	✓
Height adjustment (approximate, without mattress): Minimum: 420 mm to Maximum: 730 mm.	✓	✓
Four section uniformly perforated CRCA sheet top	✓	✓
Top and base frame made from mild steel tubes	✓	✓
Backrest adjustment	✓	✓
Height adjustment	✓	✓
Upper- leg adjustment	✓	✓
Trendelenburg & Reverse Trendelenburg adjustments	✓	✓
Auto Contour adjustment (Backrest and upper leg simultaneously adjusted)	✓	✓
User Friendly Handset	✓	✓
Attendant Control Panel (ACP) Box	✓	X
One button Cardiac Chair Position	✓	X
One button Emergency Trendelenburg	✓	X
One button Emergency Electric CPR	✓	X
Locking of certain adjustments	✓	X
Manual Emergency CPR with Quick Release Mechanism on both sides of the bed	✓	✓
Removable polymer molded head & foot boards	✓	✓
Polymer molded railings (Set of 4)	✓	✓
Pretreated and powder coated finish	✓	✓
Supplied in semi knocked down condition	✓	✓
STANDARD ACCESSORIES		
Heavy duty telescopic Infusion Pump pole with two hooks - SS	✓	✓
125 mm dia heavy duty Castors, Two with Brake	✓	✓
Buffers at four corners	✓	✓
Urine bag holder	✓	✓
Oxygen cylinder holder	✓	✓
Provision for Lifting Pole	✓	✓
Battery backup facility with built in battery charger	✓	✓
Four section mattress covered with rexine	✓	✓
OPTIONAL ACCESSORIES AT EXTRA COST		
X-ray permeable backrest with SS cassette holder (Mi-1076)	✓	✓
Lifting Pole – SS (Mi-1077)	✓	✓
Under-bed Night Lamp facility (Mi-1092)	✓	✓
Pull-out Linen Holder (Mi-1093)	✓	✓
Bed Extension facility with mattress (Mi-1094)	✓	✓
Chart holder (Mi-1082)	✓	✓
Pillow covered with rexine (Mi-1026)	✓	✓
Provision of connection for Nurse call system (Mi-1091)	✓	✓

- Accessories, colors and designs shown in the brochure illustrations may not be the part of standard product configuration.
- We reserve the right to make changes to designs, specifications and illustrations without prior notice.
- Reproduction of this brochure in any form is liable for legal action.