

BED SOLUTIONS FOWLER BEDS

SEMI MOTORIZED
MANUALLY OPERATED



Mi-9070 AX

**Motorized Fowler Bed
with Backrest and
Upper-leg on Motors**

Introduction:

Medimek is a trusted brand in the medical fraternity over a decade. Our singular and simple philosophy is to cater to the needs of all our customers with a comprehensive range of excellent quality, user friendly hospital furniture and allied products, offered in a broad spectrum of exciting shades, aesthetic patterns and creative designs and backed by prompt after sales services.

Beds are universal and essential for ensuring successful medical treatment to the patient in any hospital and medical institution. Adding option of motorized function enhances patient comfort and reduces stress of the care givers.

Ideally suited for patient care in surgical wards and deluxe rooms, **Fowler beds** offer a crucial tier of sturdy, convenient and comfortable beds in all hospitals and medical institutions. They are a preferred option while tendering home care for patients.



Mi-8046 BX

Manually Operated Fowler Bed

with Polymer head and foot boards and Collapsible railings



Mi-8047 CX

Manually Operated Fowler Bed

with SS head and foot boards with Metal panels and Collapsible railings



Mi-8047

Manually Operated Fowler Bed

with SS head and
foot boards with
Metal panels



Mi-8048

Manually Operated Fowler Bed

with SS head and foot
boards with SS tubes



Backrest



Upper
leg



Fixed SS
handle



Dual
function
detachable
handle

Options Available :

Semi Motorized Fowler Beds:

Electrical and manual adjustments:

- Only backrest electrically adjusted
- Backrest and upper leg electrically adjusted with user friendly handset
- Any non-motorized adjustment on separate screw mechanism provided with handle
- Electrical adjustments available either on actuator with manual override system OR actuator with battery back-up facility

Manually Operated Fowler Beds:

Manual adjustments:

- Backrest and upper leg section adjustments on separate screw mechanisms provided with handles

Common Features:

- Four section uniformly perforated mild steel CRCA sheet top
- Calf section adjusted on ratchet mechanism
- Top frame made from mild steel tubes
- "H" type tubular legs fitted with rubber shoes
- Four locations on bed frame for IV pole
- Pretreated and powder coated finish
- Supplied in SKD condition

Head and Foot boards :

- Polymer head and foot boards / Stainless steel head and foot boards with Metal panels / Stainless steel head and foot boards with SS tubes

Bed dimensions :

Option 1:

- Bed frame approx. dimension : 2080 mm L x 915 mm W x 600 mm H **[82"L x 36"W]**
- Overall approx. dimension (buffer-to-buffer) : 2180 mm L x 1015 mm W

Option 2:

- Bed frame approx. dimension : 1980 mm L x 840 mm W x 600 mm H **[78"L x 33"W]**
- Overall approx. dimension (buffer-to-buffer) : 2080 mm L x 940 mm W

Bed frame dimension varies with Polymer railings

Optional Accessories:

Mi-1004

Polymer railings without degree indicator (set of four)

Mi-1007

Collapsible railings with SS tubes and MS brackets (pair)

Mi-1008

Swing away type SS tubes railings (pair)

Mi-1021

Four section 100 mm thick PU foam mattress covered with rexine

Mi-1026

Pillow covered with rexine

Mi-1067

Telescopic IV pole with two hooks - SS

Mi-1068

IV rod with two hooks - SS

Mi-1059

Four heavy duty 125 mm dia castors, two with brake

Mi-6060

Crib with Attachment

Mi-1082

Chart Holder

- Accessories, colors and designs shown in the brochure illustrations may not be the part of standard product configuration. • We reserve the right to make changes to designs, specifications and illustrations without prior notice. • Reproduction of this brochure in any form is liable for legal action.