

BED SOLUTIONS RECOVERY BEDS

SEMI MOTORIZED
MANUALLY OPERATED



Mi-9040 AX
Semi Motorized
Four Section Recovery Bed
with Backrest and
Upper-leg on Motors

Introduction:

Medimek is a trusted brand in the medical fraternity over a decade. Our singular and simple philosophy is to cater to the needs of all our customers with a comprehensive range of excellent quality, user friendly hospital furniture and allied products, offered in a broad spectrum of exciting shades, aesthetic patterns and creative designs and backed by prompt after sales services.

Beds are universal and essential for ensuring successful medical treatment to the patient in any hospital and medical institution. Adding option of motorized function enhances patient comfort and reduces stress of the care givers.

Conducive to post-operative care, **Recovery Bed** is a dynamic product having the perfect functionalities for stabilizing the patient's health condition. Therefore preferred in post operative recovery rooms and emergency wards as an alternate option.



Mi-8016 BX

**Manually Operated
Four Section
Recovery Bed**
with Polymer head and
foot boards and
Collapsible railings



Mi-8017 CX

**Manually Operated
Four Section
Recovery Bed**
with SS head and foot
boards with Metal panels
and Collapsible railings



Mi-8036 BX

**Manually Operated
Two Section
Recovery Bed**
with Polymer head and
foot boards and
Collapsible railings



Mi-8037 CX

**Manually Operated
Two Section
Recovery Bed**
with SS head and foot
boards with Metal panels
and Collapsible railings



Backrest



Upper
leg



Trende-
lenburg



Reverse
Trende-
lenburg



Fixed SS
handle



Dual
function
detachable
handle

Options Available :

Semi Motorized Recovery Beds:

Electrical and manual adjustments:

- Only backrest electrically adjusted in Two Section Recovery Bed
- Only backrest electrically adjusted in Four Section Recovery Bed
- Backrest and upper leg electrically adjusted in Four Section Recovery Bed with user friendly handset
- Trendelenburg and reverse trendelenburg and any non-motorized adjustments on separate screw mechanisms provided with handles
- Electrical adjustments available either on actuator with manual override system OR actuator with battery back-up facility

Manually Operated Recovery Beds:

Manual adjustments:

- Backrest, upper leg and trendelenburg and reverse trendelenburg adjustments on separate screw mechanisms provided with handles

Common Features:

- Four / Two sections uniformly perforated mild steel CRCA sheet top
- Top frame and base frame made from mild steel tubes
- Telescopic stainless steel IV pole with four locations
- Four heavy duty castors 125 mm dia, two with brake
- Rubber buffers at four corners
- Pretreated and powder coated finish
- Supplied in SKD condition

Head and Foot boards :

- Polymer head and foot boards / SS head and foot boards with Metal panels

Four Section Recovery Beds:

Bed dimensions :

Option 1:

- Bed frame approx. dimension : 2080 mm L x 915 mm W x 600 mm H **[82" L x 36" W]**
- Overall approx. dimension (buffer-to-buffer) : 2180 mm L x 1015 mm W

Option 2:

- Bed frame approx. dimension : 1980 mm L x 840 mm W x 600 mm H **[78" L x 33" W]**
- Overall approx. dimension (buffer-to-buffer) : 2080 mm L x 940 mm W

Two Section Recovery Bed:

Bed dimensions :

Option 1:

- Bed frame approx. dimension : 1980 mm L x 915 mm W x 600 mm H **[78" L x 36" W]**
- Overall approx. dimension (buffer-to-buffer) : 2080 mm L x 1015 mm W

Option 2:

- Bed frame approx. dimension : 1905 mm L x 840 mm W x 600 mm H **[75" L x 33" W]**
- Overall approx. dimension (buffer-to-buffer) : 2005 mm L x 940 mm W

Bed frame dimension varies with polymer railings (Applies for all types of beds)

- Accessories, colors and designs shown in the brochure illustrations may not be the part of standard product configuration. • We reserve the right to make changes to designs, specifications and illustrations without prior notice. • Reproduction of this brochure in any form is liable for legal action.

Optional Accessories:

Mi-1002

Polymer railings with degree indicator (set of four)

Mi-1004

Polymer railings without degree indicator (set of four)

Mi-1007

Collapsible railings with SS tubes and MS brackets (pair)

Mi-1021

Four Section 100 mm thick PU foam mattress covered with rexine

Mi-1022

Two section 100 mm thick PU foam mattress covered with rexine

Mi-1026

Pillow covered with rexine

Mi-1066

Heavy duty Telescopic Infusion Pumps pole - SS

Mi-1077

Lifting Pole -SS

Mi-1082

Chart Holder