

BED

SOLUTIONS

SEMI FOWLER BEDS

MOTORIZED
MANUALLY OPERATED



Mi-9081 BX

**Motorized Semi Fowler Bed
with Backrest on Motor**

Introduction:

Medimek is a trusted brand in the medical fraternity over a decade. Our singular and simple philosophy is to cater to the needs of all our customers with a comprehensive range of excellent quality, user friendly hospital furniture and allied products, offered in a broad spectrum of exciting shades, aesthetic patterns and creative designs and backed by prompt after sales services.

Beds are universal and essential for ensuring successful medical treatment to the patient in any hospital and medical institution. Adding option of motorized function enhances patient comfort and reduces stress of the care givers.

Suitable in Medical wards and semi special rooms, **Semi- Fowler beds** are synonymous with the sturdiness, convenience and comfort offered to recovering patients. Also favoured by families who home care patients.



Mi-8056 BX

Manually Operated Semi Fowler Bed with Polymer head and foot boards and Collapsible railings



Mi-8057 CX

Manually Operated Semi Fowler Bed with SS head and foot boards with Metal panels and Collapsible railings



Mi-8057

Manually Operated Semi Fowler Bed with SS head and foot boards with Metal panels



Mi-8058

Manually Operated Semi Fowler Bed with SS head and foot boards with SS tubes



Backrest



Fixed SS handle



Dual function detachable handle

Options Available :

Motorized Semi Fowler Beds:

Electrical and Manual adjustment:

- Backrest electrically adjusted with user friendly handset
- Electrical adjustment available either on actuator with manual override system OR actuator with battery back-up facility

Manually Operated Semi Fowler Beds:

Manual adjustment:

- Backrest section adjustment on separate screw mechanism provided with handle

Common Features:

- Two sections uniformly perforated mild steel CRCA sheet top
- Top frame made from mild steel tubes
- "H" type tubular legs fitted with rubber shoes
- Four locations on bed frame for IV pole
- Pretreated and powder coated finish
- Supplied in SKD condition

Head and Foot boards:

- Polymer head and foot boards / Stainless steel head and foot boards with Metal panels / Stainless steel head and foot boards with SS tubes

Bed dimensions :

Option 1:

- Bed frame approx. dimension : 1980 mm L x 915 mm W x 600 mm H [78" L x 36"W]
- Overall approx. dimension (buffer-to-buffer) : 2030 mm L x 960 mm W

Option 2:

- Bed frame approx. dimension : 1905 mm L x 840 mm W x 600 mm H [75" L x 33"W]
- Overall approx. dimension (buffer-to-buffer) : 1955 mm L x 890 mm W

Bed frame dimension varies with Polymer railings

Optional Accessories:

Mi-1004

Polymer railings without degree indicator (set of four)

Mi-1007

Collapsible railings with SS tubes and MS brackets (pair)

Mi-1008

Swing away type SS tube railings (pair)

Mi-1022

Two section 100 mm thick PU foam mattress covered with rexine

Mi-1026

Pillow covered with rexine

Mi-1067

Telescopic IV pole with two hooks – SS

Mi-1068

IV rod with two hooks - SS

Mi-1059

Four heavy duty 125 mm dia castors, two with brake

Mi-6060

Crib with attachment

Mi-1082

Chart Holder

- Accessories, colors and designs shown in the brochure illustrations may not be the part of standard product configuration. • We reserve the right to make changes to designs, specifications and illustrations without prior notice. • Reproduction of this brochure in any form is liable for legal action.