

OT AND ICU SOLUTIONS



Monitor Trolleys Crash Carts

Introduction:

Medimek is a trusted brand in the medical fraternity over a decade. Our singular and simple philosophy is to cater to the needs of all our customers with a comprehensive range of excellent quality, user friendly hospital furniture and allied products, offered in a broad spectrum of exciting shades, aesthetic patterns and creative designs and backed by prompt after sales services.

One of the most essential departments of any Surgical and Medical hospital is Operation Theatre (**OT**) and Intensive Care Unit (**ICU**). Multiple clinical procedures are being performed in these departments requiring ease and swiftness of man and machines.

Monitor Trolleys serve the purpose of arranging assorted medical equipment in a useful manner and **Crash carts** are the apt trolleys for storing medicines, injections, IV fluids, report files etc. in a cogent manner.



Mi-5026



Mi-5026 A

Mi-5026 A

Monitor Trolley - SS Framework

Specification:

- Overall approx. dimension : 950 mm L x 650 mm D x 1500 mm H
- Five shelves approx. dimension: 610 mm L x 510 mm D
- Drawer approx. dimension: 580 mm L x 450 mm D
- Construction made from stainless steel tubes and sheets

Alternate:

Mi-5026

Monitor Trolley - MS Framework

- Electrical panel board having eight sockets and switches
- Oxygen Cylinder Holder - MS powder coated
- Four heavy duty castors 125 mm dia, two with brake

Optional accessory:

Mi-1090 - Scope Holder



Mi-5029

Crash Cart (Small) - MS Framework

Specification:

- Overall approx. dimension: 460 mm L x 500 mm W x 1630 mm H
- Mild steel tubular frame work and MS Shelves
- Three removable plastic bins and two polystyrene lockable storage units with two drawers having containers of different sizes
- Four heavy duty castors 100 mm dia, two with brake
- Complete with four corner buffers and SS IV Rod
- Pretreated and powder coated finish

Alternate:

Mi-5029 A

Crash Cart (Small) - SS Framework



Mi-5028 A

Crash Cart - SS Framework

Specification:

- Overall approx. dimension: 970 mm L x 500 mm W x 1630 mm H
- Stainless steel tubular framework and SS shelves
- Six removable plastic bins and two polystyrene lockable storage units with two drawers having containers of different sizes
- Four heavy duty castors 125 mm dia, two with brake
- Complete with four corner buffers, SS IV Rod, Cardiac massage board and Oxygen cylinder holder (MS)

Alternate:

Mi-5028

Crash Cart - MS Framework



Mi-5027

Wall Mounted Monitor Stand

Specification:

- Tray approx. dimension: 380 mm L x 270 mm W
- Suitably mounted on stand and fitted to the wall
- Pretreated and powder coated finish



Mi-5030

Modified Crash Cart - MS Framework

Specification:

- Overall approx. dimension:
1135 mm L x 765 mm W x 1600 mm H
- Mild steel construction of tubes and sheets
- Two rows each of six removable plastic bins in different colors
- Five MS drawers of different heights having two with multiple compartments and all drawers lockable with common lock

Alternate:

Mi-5030 A

Modified Crash Cart - SS Framework



- Four heavy duty castors 125 mm dia, two with brake
- Two SS Shelves, SS Handle with soft grip, SS Waste bin, SS Cable hanging rod and Cardiac massage board
- Four Corner buffers, two Multiuse pouches, Sliding Writing board, Reports file holder and Oxygen cylinder holder
- Available in different colour combinations
- Pretreated and powder coated finish

• Accessories, colors and designs shown in the brochure illustrations may not be the part of standard product configuration. • We reserve the right to make changes to designs, specifications and illustrations without prior notice. • Reproduction of this brochure in any form is liable for legal action.

Medimek Industries



B-50 Additional MIDC, Anand Nagar, Ambarnath [E] - 421 506
Dist. Thane, Maharashtra, India. • Tel. : +91- 251- 2621124 / 2621242 / 8888 2000 65
Email: info@medimek.com • www.medimek.com
Certifications: ISO-9001 • ISO-14001 • OHSAS-18001 • ISO-13485 • CE Certificate