

OT AND ICU SOLUTIONS



Operation Tables SS Scrubs

Introduction:

Medimek is a trusted brand in the medical fraternity over a decade. Our singular and simple philosophy is to cater to the needs of all our customers with a comprehensive range of excellent quality, user friendly hospital furniture and allied products, offered in a broad spectrum of exciting shades, aesthetic patterns and creative designs and backed by prompt after sales services.

One of the most essential departments of any Surgical and Medical hospital is Operation Theatre (**OT**) and Intensive Care Unit (**ICU**). Multiple clinical procedures are being performed in these departments requiring ease and swiftness of man and machines.

This range of **Manual Operation Tables**, provide the most comfortable and necessary surgical positioning needed during various operations. **Stainless Steel Scrubs** are also essential and indispensable products in the OT zone, for pre and post operative hand-wash for the surgical team.



Lateral Tilt



Kidney Bridge

Mi-5001

Manual Operation Table - Major (Height Adjustable)

Specification:

- Overall approx. dimension: 1830 mm L x 710 mm W
- Height adjustment on screw mechanism: 735 mm to 970 mm (approx.)
- Stainless steel sectional top and stainless steel tubular base with MS Linkages
- Removable head section, folding and removable foot section
- Lateral tilt, Kidney Bridge and Trendelenburg and reverse trendelenburg adjustments on screw mechanisms
- Complete with Padded Lithotomy Crutches (pair) with Radial Setting Clamp, Padded Arm Boards (pair), SS Removable Side Railing (single), Telescopic I V Pole, Sectional Mattress, SS Drainage Tray and Single Step Stool (MS)



Mi-5003

Manual Operation Table – Minor (Height Adjustable)

Specification:

- Overall approx. dimension: 1830 mm L x 710 mm W
- Height adjustment on screw mechanism: 735 mm to 970 mm (approx.)
- Stainless steel sectional top and stainless steel tubular base with MS Linkages
- Folding and removable foot section
- Trendelenburg and reverse trendelenburg adjustments on screw mechanism
- Complete with SS Lithotomy Rods with Stirrups (pair), Padded Arm Boards (pair), Telescopic I V Pole, SS Drainage Tray, Sectional Mattress and Single Step Stool (MS)

Mi-5015

Manual Operation Table – Ophthalmic (Height Adjustable)

Specification:

- Overall approx. dimension: 1830 mm L x 690 mm W
- Height adjustment on screw mechanism: 735 mm to 970 mm (approx.)
- Stainless steel sectional top and stainless steel tubular base with MS Linkages
- Horse-shoe shaped headrest
- Trendelenburg and reverse trendelenburg adjustments on screw mechanism
- Complete with SS Removable Side Railings, Telescopic IV Pole, Adjustable Instrument Tray, Sectional Mattress and Single Step Stool (MS)





Mi-5010

Manual Operation Table – Minor (Fixed Height)

Specification:

- Overall approx. dimension: 1830 mm L x 710 mm W x 810 mm H
- Stainless steel sectional top and stainless steel tubular base with MS Linkages
- Folding and removable foot section

- Trendelenburg and reverse trendelenburg adjustments on screw mechanism
- Complete with SS Lithotomy Rods with Stirrups (pair), Padded Arm Boards (pair), Telescopic I V Pole, SS Drainage Tray, Sectional Mattress and Single Step Stool (MS)



Foot Section Folding

Mi-5070

Wall Mounted Scrub

Specification:

- Overall approx. dimension: 480 mm W x 455 mm D
- Overall approx. Length: As specified by the customer
- Stainless steel construction with SS Washing Tray

Optional accessories



Mi-5075

Foot Operated Scrub

Specification:

- Overall approx. dimension: 555 mm W x 455 mm D
- Overall approx. Length: As specified by the customer
- Stainless steel construction with SS Washing Tray



Mi-1099

Stainless steel Slipper Stand

- Overall approx. dimension: 265 mm W x 305 mm H
- Overall approx. Length: As specified by the customer

Mi-1098

Long handle tap

- Accessories, colors and designs shown in the brochure illustrations may not be the part of standard product configuration.
- We reserve the right to make changes to designs, specifications and illustrations without prior notice.
- Reproduction of this brochure in any form is liable for legal action.

Medimek Industries



B-50 Additional MIDC, Anand Nagar, Ambarnath [E] - 421 506
Dist. Thane, Maharashtra, India. • Tel. : +91- 251- 2621124 / 2621242 / 8888 2000 65
Email: info@medimek.com • www.medimek.com
Certifications: ISO-9001 • ISO-14001 • OHSAS-18001 • ISO-13485 • CE Certificate