

OT AND ICU SOLUTIONS



**Stands
Cabinets
Trolleys**

Introduction:

Medimek is a trusted brand in the medical fraternity over a decade. Our singular and simple philosophy is to cater to the needs of all our customers with a comprehensive range of excellent quality, user friendly hospital furniture and allied products, offered in a broad spectrum of exciting shades, aesthetic patterns and creative designs and backed by prompt after sales services.

One of the most essential departments of any Surgical and Medical hospital is Operation Theatre (**OT**) and Intensive Care Unit (**ICU**). Multiple clinical procedures are being performed in these departments requiring ease and swiftness of man and machines.

Stands, Cabinets and Trolleys are utility products serving different purposes, which are practically important and assist diverse processes in these sections.



Mi-5054

Instrument Cabinet - MS Framework

Specification:

- Overall approx. dimension: 610 mm L x 355 mm D x 1535 mm H
- Mild steel tubular construction
- MS sheet panel on top, bottom and back of the cabinet
- Door, two side walls and five shelves of 5mm thick fine glass
- Lockable door with handle, pivoted to the frame
- Four legs fitted with PVC bush
- Pretreated and powder coated finish

Alternate:

Mi-5054 A

Instrument Cabinet - SS Framework

Mi-5061

Baby Resuscitation Trolley - SS Framework

Specification:

- Overall approx. dimension : 770 mm L x 415 mm W x 1030 mm H
- Stainless steel tubular construction with SS baby tray
- Complete with modular lockable storage unit with three drawers, SS Shelves, Sliding Utility board, SS IV Rod, Examination Lamp and Oxygen cylinder holder
- Four heavy duty castors 100 mm dia, two with brake





Mi-5055

Multiple Bins Trolley (Big) - MS Framework

Specification:

- Overall approx. dimension: 1370 mm L x 550 mm W x 1410 mm H
- Mild steel tubular construction
- Six rows of different size plastic bins
- Two SS trays, SS cable hanger and four corner buffers
- Four heavy duty castors 125 mm dia, two with brake
- Pretreated and powder coated finish
- Supplied in SKD condition

Alternate:

Mi-5055 A

Multiple Bins Trolley (Big) - SS Framework



Mi-5056

Multiple Bins Trolley (Small) - MS Framework

Specification:

- Overall approx. dimension: 720 mm L x 550 mm W x 1410 mm H
- Mild steel tubular construction
- Six rows of different size plastic bins
- Two SS trays, SS hanger rod and four corner buffers
- Four heavy duty castors 125 mm dia, two with brake
- Pretreated and powder coated finish
- Supplied in SKD condition

Alternate:

Mi-5056 A

Multiple Bins Trolley (Small) - SS Framework



Mi-5051

Kick Bucket - SS Framework

Specification:

- Overall approx. height: 330 mm
- Stainless steel tubular framework
- Stainless steel bowl
- Three castors 50 mm dia, without brake

Mi-5046

Soiled Linen Trolley - SS Framework

Specification:

- Overall approx. dimension: 890 mm H x 510 mm dia
- Stainless steel tubular framework
- Three heavy duty castors 100 mm dia, without brake
- Supplied with one bag of synthetic material



Mi-5052

Double Bowl Stand - SS Framework

Specification:

- Overall approx. height: 840 mm
- Stainless steel tubular framework
- Two stainless steel bowls
- Three castors 50 mm dia, without brake

Alternate:

Mi-5053

Single Bowl Stand - SS Framework



Mi-5050

Linen Serving Trolley - SS Framework

Specification:

- Overall approx. dimension:
- 1260 mm L x 600 mm W x 1120 mm H
- Stainless steel tubular construction
- Four heavy duty castors 125 mm dia, two with brake
- Complete with three SS trays for Fresh linen, SS handle, one bag of synthetic material for used linen

• Accessories, colors and designs shown in the brochure illustrations may not be the part of standard product configuration. • We reserve the right to make changes to designs, specifications and illustrations without prior notice. • Reproduction of this brochure in any form is liable for legal action.

Medimek Industries



B-50 Additional MIDC, Anand Nagar, Ambarnath [E] - 421 506
Dist. Thane, Maharashtra, India. • Tel. : +91- 251- 2621124 / 2621242 / 8888 2000 65
Email: info@medimek.com • www.medimek.com
Certifications: ISO-9001 • ISO-14001 • OHSAS-18001 • ISO-13485 • CE Certificate