

PATIENT TRANSFER SOLUTIONS



Stretcher Trolleys

Introduction:

Medimek is a trusted brand in the medical fraternity over a decade. Our singular and simple philosophy is to cater to the needs of all our customers with a comprehensive range of excellent quality, user friendly hospital furniture and allied products, offered in a broad spectrum of exciting shades, aesthetic patterns and creative designs and backed by prompt after sales services.

Patient Transfer is an indispensable and round the clock activity in every medical institution, however small or large. Smooth, safe and speedy movement is vital and of paramount importance in transferring the patient in distress or trauma. Equally, effortless transfer is crucial for care givers too.

Most suitable for patient transfer in general wards and casualty, the **Stretcher Trolley** is a strong and sturdy product. Its competence for smooth mobility and ease in patient transfer is unique. Multiple need based options are offered in this group serving different utilities.



Mi-4090



Mi-4081

Stretcher on Trolley (MS Framework and SS Stretcher Top)

Specification:

- Overall approx. dimension:
2030 mm L x 620 mm W x 880 mm H
- Stretcher top approx. dimension:
2030 mm L x 520 mm W
- Strong mild steel tubular construction

- Removable SS stretcher top supported on tubular frame
- Four heavy duty castors 150 mm dia, two with brake
- Pretreated and powder coated finish
- Supplied in SKD condition

Alternate:

Mi-4081 A

Stretcher on Trolley (SS Framework and SS Stretcher Top)



Mi-4087

Patient Shifting Trolley - One End Folding, Tiltable Top and SS Roller (MS Framework)

Specification:

- Overall approx. dimension:
1830 mm L x 600 mm W x 825 mm H
- Stretcher top approx. dimension:
1830 mm L x 500 mm W
- Strong mild steel tubular construction
- One end folding bakelite laminated board stretcher top fitted with smooth stainless steel Roller
- Stretcher top tiltable on one side
- Pair of removable stainless steel railings
- Four heavy duty castors 150 mm dia, two with brake
- Pretreated and powder coated finish

Alternate:

Mi-4088

Patient Shifting Trolley - Both Ends Folding, Tiltable Top and SS Roller (MS Framework)





Mi-4083

Stretcher Trolley - Both Ends Folding (MS Framework and SS Stretcher Top)

Specification:

- Overall approx. dimension:
1860 mm L x 600 mm W x 780 mm H (Fully extended)
- Overall approx. dimension
1000 mm L x 600 mm W x 780 mm H (Both ends folded)
- Strong mild steel tubular construction

- Fixed stainless steel stretcher top supported on tubular frame
- Both ends folding stretcher top
- Four heavy duty castors 150 mm dia, two with brake
- Pretreated and powder coated finish

Alternate:

Mi-4083 A

Stretcher Trolley - Both Ends Folding (SS Framework and SS Stretcher Top)

Optional accessories for Stretcher Trolleys:

- Mi-1080** - Oxygen cylinder holder (MS)
- Mi-1012** - Swing away type railing with stainless steel tubes (pair)
- Mi-1012 A** - Swing away type railing with mild steel tubes (pair)
- Mi-1037** - Suitable 40 mm thick PU foam mattress covered with rexine for Plain / Sectional stretcher top (**Mi-1037 A / Mi-1037 B**)
- Mi-1068** - IV rod with two hooks - SS



• Accessories, colors and designs shown in the brochure illustrations may not be the part of standard product configuration. • We reserve the right to make changes to designs, specifications and illustrations without prior notice. • Reproduction of this brochure in any form is liable for legal action.

Medimek Industries

B-50 Additional MIDC, Anand Nagar, Ambarnath [E] - 421 506

Dist. Thane, Maharashtra, India. • Tel. : +91- 251- 2621124 / 2621242 / 8888 2000 65

Email: info@medimek.com • www.medimek.com

Certifications: ISO-9001 • ISO-14001 • OHSAS-18001 • ISO-13485 • CE Certificate