

WARDCARE SOLUTIONS



Assorted Products

Introduction:

Medimek is a trusted brand in the medical fraternity over a decade. Our singular and simple philosophy is to cater to the needs of all our customers with a comprehensive range of excellent quality, user friendly hospital furniture and allied products, offered in a broad spectrum of exciting shades, aesthetic patterns and creative designs and backed by prompt after sales services.

Wardcare is an omnipresent product family in all institutions ranging from small Nursing Homes to Super speciality hospitals. **Assorted Products** are prerequisite products for handling patients and attending to their different needs comfortably.



Mi-6057

Ward Round Trolley - MS Framework

Specification:

- Overall approx. dimension: 860 mm L x 600 mm W x 1140 mm H
- Mild steel tubular framework
- Top SS Tray approx. dimension: 690 mm L x 510 mm W
- Three drawers having locks and one SS shelf
- Two tier files holder racks
- Rubber padded SS handle and four corner buffers
- Four castors 100 mm dia, two with brake
- Pretreated and powder coated finish



Mi-6058

Food Serving Trolley - SS Framework

Specification:

- Overall approx. dimension: 985 mm L x 440 mm W x 1040 mm H
- Stainless steel tubular framework
- Top SS Tray approx. dimension: 850 mm L x 400 mm W
- Double doors on both sides
- Two rows for Food tray holders
- Two trays for bottle storage
- Rubber padded SS handle
- Four castors 100 mm dia, two with brake



Mi-6074

Wheel Chair (Lifting Type) - MS Framework

Specification:

- Overall approx. dimension: 560 mm L x 890 mm W x 960 mm H
- Mild steel tubular framework with MS seat and back
- Complete with MS Foot rests, Foldable lifting handles and cushioned hand rest
- Four heavy duty castors 125 mm dia, two with brake
- Pretreated and powder coated finish

Alternate:

Mi-6073

Wheel Chair (Lifting Type) - SS Framework

Optional: Cushioned Seat and Back



Mi-6074



Mi-6073



Mi-6062



Mi-6063

Mi-6062

Locker Box Attachment for Bed

Specification:

- Overall approx. dimension: 305 mm L x 305 mm W x 295 mm H
- Mild steel Locker box with SS top and MS railings
- Box door with three louvers, handle and magnetic latch
- Pretreated and powder coated finish

Alternate:

Mi-6063

Storage Box Attachment for Bed



Mi-6060

Crib with Attachment



Mi-6061

Crib with Trolley

- Overall approx. dimension : 945 mm L x 425 mm W X 1030 mm H
- Trolley fitted with four castors 50 mm dia, two with brake

Common Specification:

- Crib approx. dimension: 720 mm L x 390 mm W x 260 mm H
- Crib made from mild steel tubes, rods and perforated sheet
- Crib swings on two vertical tubes
- Pretreated and powder coated finish

Mi-6071

Bed Side Screen (Without Curtain Cloth)

Specification:

- Overall approx. dimension: 1730 mm H x 2360 mm W (Extended)
- Mild steel tubular framework supplied with hooks
- Middle panel 1200 mm and two side panels 570 mm wide
- Six castors 50 mm dia, without brake
- Pretreated and powder coated finish

Alternate:

Mi-6071 A

Bed Side Screen (With Curtain Cloth)



Mi-6075 A

Mi-6075 B

Mi-6075 C

Mi-6075 A

Oxygen Cylinder Trolley (Pull Type)

Three castors 50 mm dia, without brake

Mi-6075 B

Oxygen Cylinder Trolley (Push Type)

Two heavy duty castors 100 mm dia, without brake

Mi-6075 C

Jumbo Oxygen Cylinder Trolley (Push Type)

Two heavy duty castors 125 mm dia, without brake

Common Specification:

- Mild steel tubular framework
- Pretreated and powder coated finish

• Accessories, colors and designs shown in the brochure illustrations may not be the part of standard product configuration. • We reserve the right to make changes to designs, specifications and illustrations without prior notice. • Reproduction of this brochure in any form is liable for legal action.

Medimek Industries



B-50 Additional MIDC, Anand Nagar, Ambarnath [E] - 421 506

Dist. Thane, Maharashtra, India. • Tel. : +91- 251- 2621124 / 2621242 / 8888 2000 65

Email: info@medimek.com • www.medimek.com

Certifications: ISO-9001 • ISO-14001 • OHSAS-18001 • ISO-13485 • CE Certificate