

WARDCARE SOLUTIONS



**Saline Stands
Stools
Trolleys**

Introduction:

Medimek is a trusted brand in the medical fraternity over a decade. Our singular and simple philosophy is to cater to the needs of all our customers with a comprehensive range of excellent quality, user friendly hospital furniture and allied products, offered in a broad spectrum of exciting shades, aesthetic patterns and creative designs and backed by prompt after sales services.

Wardcare is an omnipresent product family in all institutions ranging from small Nursing Homes to Super speciality hospitals. Saline stands, Stools and Trolleys are the quintessential, facilitating products providing fine medical care as an everyday mission.



Mi-6049

Mi-6046

Saline Stand - SS Framework

Specification:

- Height adjustment (approx.): 1350 mm to 2410 mm
- Stainless steel tubular framework
- Five pronged tubular base
- Five castors 50 mm dia, without brake
- Stainless steel rod with two hooks

Alternate:

Mi-6047

Saline Stand - MS Framework



Mi-6046

Mi-6047

Mi-6049

Syringe Pump Stand - SS Framework

Specification:

- Height adjustment (approx.): 1350 mm to 2450 mm
- Stainless steel framework with five pronged tubular base
- Complete with four adjustable SS trays and SS IV Rod
- Five 50 mm dia castors, without brake

Mi-6056

Multiple Apron Hanger - SS Framework

Specification:

- Stainless steel tubular framework
- Stainless steel hangers
- Base fitted with rubber bush

Alternate:

Mi-6055

Single Apron Hanger - SS Framework



Mi-6055



Mi-6056



Mi-6051

Mi-6051

ECG Machine Trolley - MS Framework

Specification:

- Overall approx. dimension: 360 mm L x 480 mm W x 780 mm H
- Mild steel tubular framework
- Top SS shelf approx. dimension: 360 mm L x 430 mm W
- Other two MS shelves approx. dimension: 300 mm L x 430 mm W
- SS Cable hanging rod and rubber padded SS handle
- Four castors 50 mm dia, without brake
- Pretreated and powder coated finish

Alternate:

Mi-6051 A

ECG Machine Trolley - SS Framework

Mi-6052

ECG Machine Trolley with Drawer - MS Framework



Mi-6052



Mi-6053

Mi-6053

Drug Trolley (Big) - MS Framework

Specification:

- Overall approx. dimension: 490 mm L x 580 mm W x 1000 mm H
- Mild steel tubular framework
- Top SS shelf approx. dimension: 485 mm L x 530 mm W
- Other three MS shelves approx. dimension: 485 mm L x 430 mm W
- Rubber padded SS handle
- Four castors 100 mm dia, two with brake
- Pretreated and powder coated finish

Alternate:

Mi-6053 A

Drug Trolley (Big) - SS Framework

Mi-6054

Drug Trolley (Small) - MS Framework

Overall approx. dimension:
360 mm L x 480 mm W x 780 mm H

Mi-6054 A

Drug Trolley (Small) - SS Framework

Mi-6066

Multi-purpose Stool (MS Framework and SS Top)

Specification:

- Overall approx. dimension: 335 mm L x 335 mm W X 445 mm H
- SS top approx. dimension: 280 mm L x 280 mm W
- Mild steel framework with legs fitted with PVC bush
- Pretreated and powder coated finish

Alternate:

Mi-6066 A

Multi-purpose Stool (SS Framework and SS Top)



Mi-6066

Mi-6066 A



Mi-6067

Single Step Stool - MS Framework

Specification:

- Overall approx. height: 280 mm
- Overall approx. dimension: 455 mm L x 325 mm W
- Mild steel tubular framework
- Mild steel top fitted with aluminium tread flats
- Legs fitted with rubber shoes
- Pretreated and powder coated finish

Alternate:

Mi-6067 A

Single Step Stool - SS Framework

Mi-6068

Double Step Stool - MS Framework

Specification:

- Overall approx. height: 255 mm / 490 mm
- Overall approx. dimension: 455 mm L x 595 mm W
- Mild steel tubular framework
- Mild steel top fitted with aluminium tread flats
- Legs fitted with rubber shoes
- Pretreated and powder coated finish

Alternate:

Mi-6068 A

Double Step Stool - SS Framework



- Accessories, colors and designs shown in the brochure illustrations may not be the part of standard product configuration. • We reserve the right to make changes to designs, specifications and illustrations without prior notice. • Reproduction of this brochure in any form is liable for legal action.