



Mi-4077

Emergency and Recovery Trolley (Height on Screw Mechanism and Collapsible Railings)

Specification:

- Overall approx. dimension: 1870 mm L x 730 mm W
 - Stretcher top approx. dimension: 1835 mm L x 595 mm W
 - Height adjustment on screw mechanism: 640 mm to 845 mm (approx.) (without Mattress)
 - Two section removable X-ray permeable stretcher top with sliding X-ray cassette holder
 - Backrest raised on ratchet
 - Trendelenburg and reverse trendelenburg adjustments on screw mechanism
 - Pair of collapsible stainless steel tube railings
 - Stainless steel IV pole
 - Oxygen cylinder holder and Urine bag holder
 - Buffers at four corners
 - Four heavy duty castors 125 mm dia, two with brake
 - Pretreated and powder coated finish
 - Supplied in SKD condition
- Optional Accessory:**
Two section 40 mm thick PU foam mattress covered with rexine