

Introduction :

Medimek is a trusted brand in the medical fraternity over a decade. Our singular and simple philosophy is to cater to the needs of all our customers with a comprehensive range of excellent quality, user friendly hospital furniture and allied products, offered in a broad spectrum of exciting shades, aesthetic patterns and creative designs and backed by prompt after sales services.

Patient Transfer is an indispensable and round the clock activity in every medical institution, however small or large. Smooth, safe and speedy movement is vital and of paramount importance in transferring the patient in distress or trauma. Equally, effortless transfer is crucial for care givers too.

Most suitable for patient transfer in general wards and casualty, the **Stretcher Trolley** is a strong and sturdy product. Its competence for smooth mobility and ease in patient transfer is unique. Multiple need based options are offered in this group serving different utilities.



Mi-4081

Stretcher on Trolley (MS Framework and SS Stretcher Top)

Specification:

- Overall approx. dimension: 2030 mm L x 610 mm W x 880 mm H
- Stretcher top approx. dimension: 2030 mm L x 510 mm W
- Strong mild steel tubular construction
- Removable SS stretcher top supported on tubular frame
- Four heavy duty castors 150 mm dia, two with brake
- Pretreated and powder coated finish
- Supplied in SKD condition

Alternate:

Mi-4081 A

Stretcher on Trolley (SS Framework and SS Stretcher Top)

Mi-4090

