

BED SOLUTIONS



Attendant Beds

Introduction :

Medimek is a trusted brand in the medical fraternity over a decade. Our singular and simple philosophy is to cater to the needs of all our customers with a comprehensive range of excellent quality, user friendly hospital furniture and allied products, offered in a broad spectrum of exciting shades, aesthetic patterns and creative designs and backed by prompt after sales services.

Beds are universal and essential for ensuring successful medical treatment to the patient in any hospital and medical institution. Attending a recuperating patient can be a physically exhausting experience for his family and relatives in any medical institution. **Attendant Beds** are so well designed that they offer comfort and relaxation to the patient attendants, during their hospital duties.



Mi-7021

Attendant Bed - Plain

Specification:

- Overall approx. dimension (with mattress): 1840 mm L x 610 mm W x 545 mm H
- Framework made from mild steel tubes
- Plain uniformly perforated mild steel CRCA sheet top

Optional Accessory:

Suitable mattress

- "H" type tubular legs fitted with rubber bush
- Two locations on bed frame for IV rod
- Pretreated and powder coated finish



Mi-7031

Attendant Bed - Telescopic

Specification:

- Overall approx. dimension (with mattress): 1840 mm L x 560 mm W x 440 mm H (Fully extended)
- Overall approx. dimension (with mattress): 940 mm L x 560 mm W x 510 mm H (Slided in)
- Framework made from mild steel tubes

- Mild steel CRCA sheet top with raised edges
- One section is slided under the other section
- Pretreated and powder coated finish

Optional Accessory:

Suitable mattress



Mi-7026


Attendant Bed with Panels

Specification:

- Overall approx. dimension (with mattress):
1880 mm L x 645 mm W x 750 mm H
- Top frame dimension (with mattress):
1830 mm L x 610 mm W x 545 mm H
- Construction made from mild steel tubes and sheets
- Mild steel CRCA tray type sheet top with raised edges
- Decorative panel and shelf provided at head end
- Pretreated and powder coated finish

Optional Accessory:

Suitable mattress



Mi-7036

Attendant Bed with Two Cabinets

Specification:

- Overall approx. dimension (with mattress):
1840 mm L x 610 mm W x 560 mm H
- Construction made from mild steel tubes and sheets
- Mild steel CRCA tray type sheet top with raised edges
- Two cabinets with locks and open shelf in middle
- Pretreated and powder coated finish

Optional Accessory:

Suitable mattress

Mi-7046

Low Back Attendant Bed with Mattress (Backrest Reclining on Gas spring and manually operated legrest)

- Overall approx. dimension (with mattress and hand rest): 1715 mm L x 705 mm W
- Sleeping level height (with mattress): 430 mm H
- Seating, Reclining and Sleeping positions manually adjusted



Reclining position



Sleeping position

Common Specification:

- Construction made from mild steel tubes and sheets
- Supplied with suitable mattress
- Cushioned handrest
- Tubular legs fitted with rubber shoes
- Pretreated and powder coated finish

Mi-7052

High Back Attendant Bed with Mattress (Backrest Reclining on Gas spring and manually operated legrest)

- Overall approx. dimension (with mattress and hand rest): 1815 mm L x 760 mm W
- Sleeping level height (with mattress) : 430 mm H
- Seating, Reclining and Sleeping positions obtained by adjusting backrest on gas spring mechanism and manually operated leg rest
- Provision to operate gas spring from both sides



Seating position



Sleeping position



Reclining position

- Accessories, colors and designs shown in the brochure illustrations may not be the part of standard product configuration.
- We reserve the right to make changes to designs, specifications and illustrations without prior notice.
- Reproduction of this brochure in any form is liable for legal action.