

GYNAEC AND OBSTETRIC SOLUTIONS



Motorized Birthing Bed
Motorized Gynaecology Chair
Obstetric Labour Table

Introduction :

Medimek is a trusted brand in the medical fraternity over a decade. Our singular and simple philosophy is to cater to the needs of all our customers with a comprehensive range of excellent quality, user friendly hospital furniture and allied products, offered in a broad spectrum of exciting shades, aesthetic patterns and creative designs and backed by prompt after sales services.

In this family, the ergonomically designed **Motorized Gynaecology Chair** offers ease of access for patient examination and procedures. **Motorized Birthing Bed** is designed to deliver a soothing and dignified birthing experience to mothers and efficient functionality to nursing staff, in addition to enhanced mobility. The popular option of stainless steel **Obstetric Labour Table** offers convenient and comfortable positioning and utility for safe child birth.



Mi-7075

Motorized Gynaecology Chair

Specification:

- Overall approx. dimension: 1675 mm L x 660 mm W
- Height adjustment (approx.): 630 mm to 940 mm
- Backrest, height, pelvic tilt adjustments electrically operated on actuators with user friendly handset
- Leg rest adjustment manually operated
- Complete with SS tray, Padded Lithotomy Crutches (pair)
- PU foam upholstered top with perennial "U" cut covered with rexine
- Six PVC Levellers
- Pretreated and powder coated finish
- Supplied in SKD condition

Optional Accessory:

Battery back-up



Mi-7077

Motorized Birthing Bed

Specification:

- Overall approx. dimension: 2010 mm L x 1020 mm W
- Height adjustment (approx. without mattress): 535 mm to 900 mm
- Three section laminated top with perennial "U" notch
- Backrest, height and trendelenburg and reverse trendelenburg adjustments electrically operated on actuators with user friendly handset
- Leg section can be telescoped under the backrest section for lithotomic positions
- Removable Polymer head and foot boards
- Pair of Polymer railings
- Middle and Leg section tops are removable for cleaning
- Complete with SS tray, SS telescopic IV Pole, SS hand grips, Padded Lithotomy Crutches (pair) and rubber buffers
- PU foam mattress covered with rexine
- Four heavy duty castors of 150 mm dia, with central and directional locking system
- Battery back-up with built-in charger
- Pretreated and powder coated finish
- Supplied in SKD condition



Mi-7085

Obstetric Labour Table – Telescopic (SS Framework)

Specification:

- Overall approx. dimension: 1930 mm L x 890 mm W x 790 mm H (Fully extended without mattress)
- Strong stainless steel tubular framework
- Three section stainless steel table top
- Backrest raised on ratchet
- Middle section with perennial "U" cut
- Leg section can be telescoped under the main frame
- Trendelenburg and reverse trendelenburg adjustments on screw with SS handle
- Removable SS head and side railings
- Complete with SS Bowl, SS Telescopic IV Rod, SS Hand grips, Padded Lithotomy Crutches (pair) and Single Step Stool (MS)
- Foam mattress covered with rexine
- Supplied in SKD condition



• Accessories, colors and designs shown in the brochure illustrations may not be the part of standard product configuration. • We reserve the right to make changes to designs, specifications and illustrations without prior notice. • Reproduction of this brochure in any form is liable for legal action.

Medimek Industries

B/50-53, Additional MIDC, Anand Nagar, Ambernath [E] - 421 506

Dist. Thane, Maharashtra, India. • Tel. : +91- 251- 2621124 / 2621242 / 8888 2000 65

Email: sales@medimek.com • www.medimek.com

Certifications: ISO-9001 • ISO-14001 • ISO-45001 • ISO-50001 • ISO-13485 • CE Certificate

