



Mi-5028 A

Crash Cart - SS Framework

Specification:

- Overall approx. dimension: 970 mm L x 480 mm W x 1655 mm H
- Stainless steel tubular framework and SS shelves
- Six removable plastic bins and two polystyrene lockable storage units with two drawers having containers of different sizes
- Four heavy duty castors 125 mm dia, two with brake
- Complete with four corner buffers, SS IV Rod, Cardiac massage board and Oxygen cylinder holder (MS)

Alternate:

Mi-5028

Crash Cart - MS Framework

Mi-5029

Crash Cart (Small) - MS Framework

Specification:

- Overall approx. dimension: 660 mm L x 480 mm W x 1630 mm H
- Mild steel tubular frame work and MS Shelves
- Three removable plastic bins and two polystyrene lockable storage units with two drawers having containers of different sizes
- Four heavy duty castors 100 mm dia, two with brake
- Complete with four corner buffers , SS IV Rod and Oxygen cylinder holder (MS)
- Pretreated and powder coated finish

Alternate:

Mi-5029 A

Crash Cart (Small) - SS Framework



Mi-5087

Wall Mounted Monitor Stand - MS

Specification:

- Tray approx. dimension: 380 mm L x 270 mm W
- Suitably mounted on stand and fitted to the wall
- PVC coated wire mesh basket
- Pretreated and powder coated finish

Mi-5030

Crash Cart (Big) - MS Framework and MS Drawers

Specification:

- Overall approx. dimension: 1040 mm L x 590 mm W x 1600 mm H
- Mild steel construction of tubes and sheets
- Two rows each of six removable plastic bins in different colors
- Five MS drawers of different heights having two with multiple compartments and all drawers lockable with common lock
- Four heavy duty castors 125 mm dia, two with brake
- Two SS Shelves, SS Handle with soft grip, SS Waste bin, SS Cable hanging rod
- Four Corner buffers, two Multiuse pouches, Sliding Writing board, Reports file holder and Oxygen cylinder holder
- Available in different colour combinations
- Pretreated and powder coated finish

Alternate:

Mi-5030 A

Crash Cart (Big) - SS Framework and SS Drawers



Mi-5024

Crash Cart (Small) - MS Framework and MS Drawers

Specification:

- Overall approx. dimension: 740 mm L x 590 mm W x 1450 mm H
- Mild steel construction of tubes and sheets
- One row of three removable plastic bins
- Four MS drawers of different heights having one with multiple compartments
- Four heavy duty castors 100 mm dia, two with brake
- Two SS Shelves, SS Handle with soft grip, Four corner buffers, Sliding MS writing board, SS IV rod, Oxygen cylinder holder (MS), SS Waste bin and Reports file holder
- Available in different colour combinations
- Pretreated and powder coated finish

Alternate:

Mi-5024 A

Crash Cart (Small) - SS Framework and SS Drawers

- Accessories, colors and designs shown in the brochure illustrations may not be the part of standard product configuration.
- We reserve the right to make changes to designs, specifications and illustrations without prior notice.
- Reproduction of this brochure in any form is liable for legal action.