

BED
SOLUTIONS
ELECTRIC ICU BEDS



Mi-9005 A
Electric ICU Bed
with Patient Handset and
Central Locking Castors

Introduction :

Medimek is a trusted brand in the medical fraternity over a decade. Our singular and simple philosophy is to cater to the needs of all our customers with a comprehensive range of excellent quality, user friendly hospital furniture and allied products, offered in a broad spectrum of exciting shades, aesthetic patterns and creative designs and backed by prompt after sales services.

Beds are universal and essential for ensuring successful medical treatment to the patient in any hospital and medical institution. Adding option of motorized function enhances patient comfort and reduces stress of the care givers.

Electric ICU Bed with Patient Handset is devised as an asset in patient care that performs effortlessly when you need maximum comfort. Along with the ease of electrically operated functional positions for the patient, it also complements the caregiver with sublime efficiency. Our endeavour is to provide high-end, user compatible products to administer advanced and multiple treatments to the patient.

Mi-9005 A

Electric ICU Bed
With Patient handset
and Central Locking Castors



Technical Specifications for Mi-9005 A & Mi-9005 B:

- Overall approximate dimension: 2200 mm L x 1060 mm W
- Mattress approximate dimension: 1960 mm L x 860 mm W
- Height adjustable (approximate, without mattress): Minimum: 430 mm to Maximum: 730 mm
- Backrest, upper leg, height adjustment, trendelenburg and reverse trendelenburg adjustments electrically operated on four electro mechanical actuators with user friendly handset

Adjustments obtained Electrically with Patient Handset:

- Back section
- Upper leg section
- Height
- Trendelenburg & Reverse Trendelenburg
- Electric CPR
- Under Bed Light (If opted for initially)

Standard Features:

- Manual Emergency CPR with instant lowering of backrest obtained with the help of quick release mechanism lever provided on both sides of the bed
- Four section uniformly perforated mild steel CRCA sheet top
- Top and base frame made from mild steel tubes
- Easily removable Polymer molded head and foot boards
- Polymer railings with degree indicator - set of four
- Buffers at four corners
- Urine bag holder
- Provision for Lifting pole
- Pretreated and epoxy powder coated finish
- Supplied in semi knocked down condition

Bed Models:

- **Mi-9005 A** : Bed supplied with 150 mm dia heavy duty central locking castors with directional lock having brake mechanism lever
- **Mi-9005 B** : Bed supplied with 125 mm dia heavy duty castors, two with brake two without brake

Optional Accessory :

Adjustments obtained Electrically with Atendant Control Panel (ACP) :

- Back section
- Upper leg section
- Height
- Trendelenburg & Reverse Trendelenburg
- Auto Contour (Backrest and upper leg simultaneously adjusted)
- One button Cardiac Chair Position
- One button Emergency Trendelenburg
- One button Emergency Electric CPR
- Under Bed Light (If opted for initially)
- Locking all or any desired adjustments
- Battery Charge Indicator

Bed variant with Bed Extension facility

Bed variant with X – ray Permeable backrest

- Accessories, colors and designs shown in the brochure illustrations may not be the part of standard product configuration. • We reserve the right to make changes to designs, specifications and illustrations without prior notice. • Reproduction of this brochure in any form is liable for legal action.

Optional Accessories:

Mi-1021

Four Section 100 mm thick PU Foam Mattress covered with rexine

Mi-1026

Pillow covered with rexine

Mi-1066

Heavy duty Telescopic Infusion Pump Pole - SS

Mi-1076

X-ray permeable backrest With cassette holder

Mi-1077

Lifting Pole – SS

Mi-1053

Oxygen Cylinder Cage Attachment for Bed

Mi-1092

Under bed Light facility

Mi-1083 B

PVC Coated Wire mesh file basket

Mi-1093

Pull-out Linen holder

Mi-1094

Bed Extension Facility with Mattress

Mi-1083 C

Urine Pot Basket

Battery Backup Facility