

BED SOLUTIONS

FOWLER BEDS

MOTORIZED

MANUALLY OPERATED



Mi-9070

**Motorized Fowler Bed
with Backrest and
Upper-leg on Motors,
Polymer head and foot boards**

Optional Accessory:

Polymer railings

Four section mattress

Castors

Introduction :

Medimek is a trusted brand in the medical fraternity over a decade. Our singular and simple philosophy is to cater to the needs of all our customers with a comprehensive range of excellent quality, user friendly hospital furniture and allied products, offered in a broad spectrum of exciting shades, aesthetic patterns and creative designs and backed by prompt after sales services.

Beds are universal and essential for ensuring successful medical treatment to the patient in any hospital and medical institution. Adding option of motorized function enhances patient comfort and reduces stress of the care givers. Ideally suited for patient care in surgical wards and deluxe rooms, **Fowler Beds** offer a crucial tier of sturdy, convenient and comfortable beds in all hospitals and medical institutions. They are a preferred option while tendering home care for patients.



Mi-8046

Manually Operated Fowler Bed

with Polymer head and foot boards

Optional Accessory:
SS Collapsible railings
Castors



Mi-8047

Manually Operated Fowler Bed

with SS head and foot boards
with Metal panels

Optional Accessory:
MS Collapsible railings
Castors



Mi-8047

Manually Operated Fowler Bed

with SS head and foot boards with Metal panels



Mi-8048

Manually Operated Fowler Bed

with SS head and foot boards with SS tubes



Backrest



Upper leg



Fixed SS handle



Detachable handle for Dual Function Acuator

Options Available :

Motorized Fowler Beds :

Electrical and manual adjustments :

- Only backrest electrically adjusted
- Backrest and upper leg electrically adjusted with user friendly handset
- Any non-motorized adjustment on separate screw mechanism provided with handle
- Electrical adjustments available on actuator with or without manual override system

Manually Operated Fowler Beds :

Manual adjustments :

- Backrest and upper leg section adjustments on separate screw mechanisms provided with handles

Common Features :

- Four section uniformly perforated mild steel CRCA sheet top
- Calf section adjusted on ratchet mechanism
- Top frame made from mild steel tubes
- "H" type tubular legs fitted with rubber shoes
- Four locations on bed frame for IV pole, Urine bag holder
- Pretreated and powder coated finish
- Supplied in SKD condition

Head and Foot boards :

- Polymer head and foot boards / Stainless steel head and foot boards with Metal panels / Stainless steel head and foot boards with SS tubes

Railings :

- Polymer railings / stainless steel collapsible railings / mild steel collapsible railings / Aluminium collapsible railings

Bed dimensions :

Option 1 :

- Bed frame approx. dimension : 2080 mm L x 915 mm W x 600 mm H **[82"L x 36"W]**
- Overall approx. dimension (buffer-to-buffer) : 2170 mm L x 1005 mm W

Option 2 :

- Bed approx. dimension : 1980 mm L x 840 mm W x 600 mm H **[78"L x 33"W]**
- Overall approx. dimension (buffer-to-buffer) : 2065 mm L x 920 mm W

Bed dimension varies with Polymer railings

Optional Accessories:

Mi-1002

Polymer railings with degree indicator - set of two

Mi-1003

Polymer railings with degree indicator - set of four

Mi-1007

Collapsible railings with SS tubes and MS brackets (pair)

Mi-1007 A

Collapsible railings with MS tubes and MS brackets (pair)

Mi-1010

Collapsible railings with Aluminum tubes (pair)

Mi-1021

Four section PU Foam Mattress covered with rexine

Mi-1026

Pillow covered with rexine

Mi-1053

Oxygen Cylinder Cage

Mi-1059

Heavy duty 125 mm dia castors, two with brake & two without brake

Mi-1066

Heavy duty Telescopic Infusion pump pole – SS

Mi-1067

Telescopic IV pole with Two Hooks – SS

Mi-1068

IV rod with Two Hooks – SS

Mi-1083

PVC Coated Wire mesh file basket

Mi-1083 C

Urine Pot Basket

- Accessories, colors and designs shown in the brochure illustrations may not be the part of standard product configuration. • We reserve the right to make changes to designs, specifications and illustrations without prior notice. • Reproduction of this brochure in any form is liable for legal action.

Medimek Industries



B/50-53, Additional MIDC, Anand Nagar, Ambarnath [E] - 421 506

Dist. Thane, Maharashtra, India. • Tel. : +91- 251- 2621124 / 2621242 / 8888 2000 65

Email: sales@medimek.com • www.medimek.com

Certifications: ISO-9001 • ISO-14001 • ISO-45001 • ISO-50001 • ISO-13485 • CE Certificate