

EXAMINATION SOLUTIONS



Gynaec and General Purpose Examination Couches Tables

Introduction :

Medimek is a trusted brand in the medical fraternity over a decade. Our singular and simple philosophy is to cater to the needs of all our customers with a comprehensive range of excellent quality, user friendly hospital furniture and allied products, offered in a broad spectrum of exciting shades, aesthetic patterns and creative designs and backed by prompt after sales services.

This indispensable and versatile product family is intrinsic to both Gynaec and General purpose patient examination rooms. **Examination Solutions** offer a comprehensive range of convenient, comfortable and durable Examination Couches and Tables that are easy to maintain, available in ergonomic designs and spectacular colour combinations with ample storage facility.



Mi-4011

Examination Couch with Head rest on Gas Spring

Specification:

- Overall approx. dimension: 1855 mm L x 620 mm W x 875 mm H
- Table top approx. dimension: 1855 mm L x 570 mm W
- Construction made from mild steel CRCA tubes and sheets
- Two section mild steel CRCA sheet top
- Headrest adjustable on gas spring
- Three drawers with locks
- Three cabinets with two having locks
- Complete with B.P. apparatus tray, Foam mattress covered with rexine and Built-in Step Stool
- Pretreated and powder coated finish

Mi-4012

Examination Couch - Plain Top with a Drawer and Four Cabinets

Specification:

- Overall approx. dimension: 1830 mm L x 580 mm W x 875 mm H
- Table top approx dimension 1830 mm L x 570 mm W
- Construction made from mild steel CRCA tubes and sheets
- Plain mild steel CRCA sheet top
- One of four cabinets and a drawer having locks
- Foam mattress covered with rexine
- Built-in Sliding Step Stool and Sliding Utility Board
- Pretreated and powder coated finish



Mi-4016

Examination Couch - Plain Top with Three Drawers and Three Cabinets

Specification:

- Overall approx. dimension: 1830 mm L x 570 mm W x 850 mm H
- Table top approx. dimension: 1830 mm L x 570 mm W
- Construction made from mild steel CRCA tubes and sheets
- Plain mild steel CRCA sheet top
- All three cabinets and three drawers having locks
- Complete with Foam mattress covered with rexine and Built-in Sliding Step Stool and Sliding Utility Board
- Pretreated and powder coated finish



Mi-4026

Gynaec Examination Couch

Specification:

- Overall approx. dimension: 1830 mm L x 580 mm W x 875 mm H
- Table top approx dimension: 1830 mm L x 570 mm W
- Construction made from mild steel CRCA tubes and sheets
- Plain mild steel CRCA sheet top with "U" cut
- One of three cabinets and a drawer having locks
- Foam mattress covered with rexine
- Built-in Sliding Step Stool and Sliding Utility Board
- Pretreated and powder coated finish



Mi-4015

Examination Couch - Plain Top with Decorative Columns

Specification:

- Overall approx. dimension: 1830 mm L x 580 mm W x 875 mm H
- Table top approx dimension: 1830 mm L x 570 mm W
- Mild steel tubular frame work
- Drawer and cabinet having locks
- Plain mild steel CRCA sheet top
- Foam mattress covered with rexine
- Built-in Sliding Utility Board and Single Step Stool
- Legs fitted with rubber shoes
- Pretreated and powder coated finish



Mi-4013

Examination Couch - Plain Top with a Drawer and a Cabinet

Specification:

- Overall approx. dimension: 1830 mm L x 580 mm W x 875 mm H
- Mild steel tubular frame work
- Drawer and cabinet having locks
- Plain mild steel CRCA sheet top
- Foam mattress covered with rexine
- Built-in Sliding Utility Board and Single Step Stool
- Legs fitted with rubber shoes
- Pretreated and powder coated finish



Mi-4014

Examination Table with Mattress

Specification:

- Overall approx. dimension: 1840 mm L x 605 mm W x 810 mm H
- Table top approximate dimension: 1840 mm L x 605 mm W
- Mild steel tubular frame work
- Two section mild steel CRCA sheet top
- Headrest raised on ratchet
- Foam mattress covered with rexine
- Single Step Stool
- Legs fitted with rubber shoes
- Pretreated and powder coated finish



Mi-4027

Gynaec Examination Table with "U" cut, Mattress, SS Lithotomy Rods and SS Bowl

Specification:

- Overall approx. dimension: 1840 mm L x 740 mm W x 810 mm H
(Other Specification as per Product Code Mi-4014)



• Accessories, colors and designs shown in the brochure illustrations may not be the part of standard product configuration. • We reserve the right to make changes to designs, specifications and illustrations without prior notice. • Reproduction of this brochure in any form is liable for legal action.

Medimek Industries

B/50-53, Additional MIDC, Anand Nagar, Ambarnath [E] - 421 506

Dist. Thane, Maharashtra, India. • Tel. : +91- 251- 2621124 / 2621242 / 8888 2000 65

Email: sales@medimek.com • www.medimek.com

Certifications: ISO-9001 • ISO-14001 • ISO-45001 • ISO-50001 • ISO-13485 • CE Certificate