

Introduction :

Medimek is a trusted brand in the medical fraternity over a decade. Our singular and simple philosophy is to cater to the needs of all our customers with a comprehensive range of excellent quality, user friendly hospital furniture and allied products, offered in a broad spectrum of exciting shades, aesthetic patterns and creative designs and backed by prompt after sales services.

Beds are universal and essential for ensuring successful medical treatment to the patient in any hospital and medical institution. Adding option of motorized function enhances patient comfort and reduces stress of the care givers.

ICU Bed is a holistic product designed with excellent features for patient positioning together with mobility, comfort and luxury. A composite hospital bed, it is acclaimed for Emergency care and Critical care because of its ability to offer all primarily important positions.



Mi-9020 B

Semi Motorized ICU Bed with Backrest and Height on Motors

with Polymer head and foot boards, Collapsible railings and Castors

Optional Accessory:

Four section mattress
Telescopic IV Pole with two hooks - SS



Mi-8001 A

Manually Operated ICU Bed

with Polymer head and foot boards, Polymer railings and Castors

Optional Accessory:

Four section mattress
Heavy duty Telescopic Infusion pump pole - SS



Mi-8001 B

Manually Operated ICU Bed
with Polymer head and foot boards,
Collapsible railings
and Castors

Optional Accessory:
Four section mattress
Telescopic IV Pole with
two hooks - SS



Mi-8002 B

Manually Operated ICU Bed
with SS head and foot boards
with Metal panels,
Collapsible railings
and Castors

Optional Accessory:
Four section mattress
Telescopic IV Pole with
two hooks - SS



Backrest



Upper
leg



Height



Trende-
lenburg



Reverse
Trende-
lenburg



Fixed SS
handle



Detachable
handle for
Dual Function
Actuator

Options Available :

Semi Motorized ICU Beds :

Electrical and Manual adjustments:

- Only backrest electrically adjusted
- Backrest and height / upper leg electrically adjusted OR
- Backrest, upper leg and height electrically adjusted with user friendly handset
- Trendelenburg and reverse trendelenburg and any non-motorized adjustments on separate screw mechanisms provided with handles
- Electrical adjustments available on actuator with or without manual override system

Manually Operated ICU Beds :

Manual adjustments:

- Backrest, upper leg, height and trendelenburg and reverse trendelenburg adjustments on separate screw mechanisms provided with handles

Height adjustment: (With Polymer Railings)

- Height adjustment (approx. without mattress): 545 mm to 740 mm

Height adjustment: (With Collapsible Railings)

- Height adjustment (approx. without mattress): 515 mm to 700 mm

Common Features :

- Four section uniformly perforated mild steel CRCA sheet top
- Top frame and base frame made from mild steel tubes
- Telescopic stainless steel IV pole with four locations
- Four heavy duty castors 125 mm dia, two with brake
- Rubber buffers at four corners, Urine bag holder
- Pretreated and powder coated finish
- Supplied in SKD condition

Head and Foot boards :

- Polymer head and foot boards / Stainless steel head and foot boards with Metal panels / Stainless steel head and foot boards with SS tubes

Railings :

- Polymer railings / stainless steel collapsible railings/mild steel collapsible railings / Aluminium collapsible railings

Bed dimensions :

Option 1 :

- Bed frame approx. dimension: 2080 mm L x 915 mm W [82"L x 36"W]
- Overall approx. dimension (buffer-to-buffer) : 2180 mm L x 1005 mm W

Option 2 :

- Bed frame approx. dimension: 1980 mm L x 840 mm W [78"L x 33"W]
- Overall approx. dimension (buffer-to-buffer) : 2080 mm L x 920 mm W

Bed dimension varies with Polymer railings

Bed variant with X- ray Permeable backrest

Optional Accessories:

Mi-1021

Four section PU Foam
Mattress covered with rexine

Mi-1026

Pillow covered with rexine

Mi-1053

Oxygen Cylinder Cage

Mi-1066

Heavy duty Telescopic
Infusion pump pole – SS

Mi-1067

Telescopic IV pole with Two Hooks – SS

Mi-1077

Lifting Pole – SS

Mi-1083

PVC Coated
Wire mesh file basket

Mi-1093

Pull-out Linen Holder

Mi-1076

X-ray permeable backrest with
SS Cassette holder

Mi-1083 C

Urine Pot Basket

- Accessories, colors and designs shown in the brochure illustrations may not be the part of standard product configuration. • We reserve the right to make changes to designs, specifications and illustrations without prior notice. • Reproduction of this brochure in any form is liable for legal action.