

## Introduction :

**Medimek** is a trusted brand in the medical fraternity over a decade. Our singular and simple philosophy is to cater to the needs of all our customers with a comprehensive range of excellent quality, user friendly hospital furniture and allied products, offered in a broad spectrum of exciting shades, aesthetic patterns and creative designs and backed by prompt after sales services.

One of the most essential departments of any Surgical and Medical hospital is Operation Theatre (**OT**) and Intensive Care Unit (**ICU**). Multiple clinical procedures are being performed in these departments requiring ease and swiftness of man and machines.

**Monitor Trolleys** serve the purpose of arranging assorted medical equipment in a useful manner and **Crash Carts** are the apt trolleys for storing medicines, injections, IV fluids, report files etc. in a cogent manner.



Mi-5026 A

### Mi-5026 A

#### Monitor Trolley - SS Framework

##### Specification:

- Overall approx. dimension :  
925 mm L x 640 mm D x 1480 mm H
- Five shelves approx. dimension: 605 mm L x 510 mm D
- Drawer approx. dimension: 580 mm L x 440 mm D
- Construction made from stainless steel tubes and sheets

- Electrical panel board having eight sockets and switches
- Oxygen Cylinder Holder - Stainless steel
- Four heavy duty castors 125 mm dia, two with brake

##### Optional accessory:

**Mi-1090** - Scope Holder

#### Alternate:

### Mi-5026

#### Monitor Trolley - MS Framework