

Introduction :

Medimek is a trusted brand in the medical fraternity over a decade. Our singular and simple philosophy is to cater to the needs of all our customers with a comprehensive range of excellent quality, user friendly hospital furniture and allied products, offered in a broad spectrum of exciting shades, aesthetic patterns and creative designs and backed by prompt after sales services.

In this family, the ergonomically designed **Motorized Gynaecology Chair** offers ease of access for patient examination and procedures. **Motorized Birthing Bed** is designed to deliver a soothing and dignified birthing experience to mothers and efficient functionality to nursing staff, in addition to enhanced mobility. The popular option of stainless steel **Obstetric Labour Table** offers convenient and comfortable positioning and utility for safe child birth.



Mi-7075

Motorized Gynaecology Chair

Specification:

- Overall approx. dimension: 1675 mm L x 660 mm W
- Height adjustment (approx.): 630 mm to 940 mm
- Backrest, height, pelvic tilt adjustments electrically operated on actuators with user friendly handset
- Leg rest adjustment manually operated
- Complete with SS tray, Padded Lithotomy Crutches (pair)
- PU foam upholstered top with perennial "U" cut covered with rexine
- Six PVC Levellers
- Pretreated and powder coated finish
- Supplied in SKD condition

Optional Accessory:

Battery back-up