


Mi-5070 A

OT Scrub - Wall Mounted

Specification:

- Overall approx. dimension: 510 mm W x 455 mm D x 990 mm H
- Overall approx. Length: As specified by the customer
- Stainless steel construction with SS Washing Try, Soap Holder and Soap Dispenser.



Mi-5075 A

OT Scrub - Foot Operated

Specification:

- Overall approx. dimension: 595 mm W x 455 mm D x 1245 mm H
- Overall approx. Length: As specified by the customer
- Stainless steel construction with SS Washing Try, Soap Holder and Soap Dispenser.

Optional accessories :

Mi-1099

OT Stand - SS

- Overall approx. dimension: 265 mm W x 310 mm
- Overall approx. Length: As specified by the customer



Mi-5077 A

OT Scrub - Stand Model

Specification:

- Overall approx. dimension: 505 mm W x 455 mm D x 1420 mm H
- Overall approx. Length: As specified by the customer
- Stainless steel construction with SS Washing Try, Soap Holder and Soap Dispenser.

• Accessories, colors and designs shown in the brochure illustrations may not be the part of standard product configuration. • We reserve the right to make changes to designs, specifications and illustrations without prior notice. • Reproduction of this brochure in any form is liable for legal action.

Medimek Industries



B/50-53, Additional MIDC, Anand Nagar, Ambarnath [E] - 421 506

Dist. Thane, Maharashtra, India. • Tel. : +91- 251- 2621124 / 2621242 / 8888 2000 65

Email: sales@medimek.com • www.medimek.com

Certifications: ISO-9001 • ISO-14001 • ISO-45001 • ISO-50001 • ISO-13485 • CE Certificate