

BED SOLUTIONS SEMI FOWLER BEDS

MOTORIZED
MANUALLY OPERATED



Mi-9081

**Motorized Semi Fowler Bed
with Backrest on Motor,
with Polymer head and foot
boards**

Optional Accessory:

Two section mattress
SS Collapsible railings
Castors

Introduction :

Medimek is a trusted brand in the medical fraternity over a decade. Our singular and simple philosophy is to cater to the needs of all our customers with a comprehensive range of excellent quality, user friendly hospital furniture and allied products, offered in a broad spectrum of exciting shades, aesthetic patterns and creative designs and backed by prompt after sales services.

Beds are universal and essential for ensuring successful medical treatment to the patient in any hospital and medical institution. Adding option of motorized function enhances patient comfort and reduces stress of the care givers.

Suitable in Medical wards and semi special rooms, **Semi- Fowler beds** are synonymous with the sturdiness, convenience and comfort offered to recovering patients. Also favoured by families who home care patients.



Mi-8056

Manually Operated Semi Fowler Bed

with Polymer head and foot boards

Optional Accessory:
SS Collapsible railings



Mi-8057

Manually Operated Semi Fowler Bed

with SS head and foot boards with Metal panels

Optional Accessory:
MS Collapsible railings



Mi-8057

Manually Operated Semi Fowler Bed

with SS head and foot boards with Metal panels



Mi-8058

Manually Operated Semi Fowler Bed

with SS head and foot boards with SS tubes



Backrest



Fixed SS handle



Detachable handle for Dual function Actuator

Options Available :

Motorized Semi Fowler Beds :

Electrical and Manual adjustment :

- Backrest electrically adjusted with user friendly handset
- Electrical adjustment available on actuator with or without manual override system

Manually Operated Semi Fowler Beds :

Manual adjustment :

- Backrest section adjustment on separate screw mechanism provided with handle

Common Features :

- Two sections uniformly perforated mild steel CRCA sheet top
- Top frame made from mild steel tubes
- "H" type tubular legs fitted with rubber shoes
- Four locations on bed frame for IV pole, Urine bag holder
- Pretreated and powder coated finish
- Supplied in SKD condition

Head and Foot boards :

- Polymer head and foot boards / Stainless steel head and foot boards with Metal panels / Stainless steel head and foot boards with SS tubes

Railings :

- Polymer railings / stainless steel collapsible railings / mild steel collapsible railings / Aluminium collapsible railings

Bed dimensions :

Option 1 :

- Bed frame approx. dimension : 1980 mm L x 915 mm W x 600 mm H [78"L x 36"W]
- Overall approx. dimension (buffer-to-buffer) : 2030 mm L x 960 mm W

Option 2 :

- Bed frame approx. dimension : 1905 mm L x 840 mm W x 600 mm H [75" L x 33"W]
- Overall approx. dimension (buffer-to-buffer) : 1955 mm L x 890 mm W

Bed dimension varies with Polymer railings

Optional Accessories:

Mi-1002

Polymer railings with degree indicator - set of two

Mi-1003

Polymer railings with degree indicator - set of four

Mi-1007

Collapsible railings with SS tubes and MS brackets (pair)

Mi-1007 A

Collapsible railings with MS tubes and MS brackets (pair)

Mi-1010

Collapsible railings with Aluminum tubes (pair)

Mi-1022

Two section PU Foam Mattress covered with rexine

Mi-1026

Pillow covered with rexine

Mi-1053

Oxygen Cylinder Cage

Mi-1059

Heavy duty 125 mm dia castors, two with brake & two without brake

Mi-1067

Telescopic IV pole with Two Hooks – SS

Mi-1068

IV rod with Two Hooks – SS

Mi-1083

PVC Coated Wire mesh file basket

Mi-1083 C

Urine Pot Basket

- Accessories, colors and designs shown in the brochure illustrations may not be the part of standard product configuration. • We reserve the right to make changes to designs, specifications and illustrations without prior notice. • Reproduction of this brochure in any form is liable for legal action.

Medimek Industries

B/50-53, Additional MIDC, Anand Nagar, Ambarnath [E] - 421 506

Dist. Thane, Maharashtra, India. • Tel. : +91- 251- 2621124 / 2621242 / 8888 2000 65

Email: sales@medimek.com • www.medimek.com

Certifications: ISO-9001 • ISO-14001 • ISO-45001 • ISO-50001 • ISO-13485 • CE Certificate

