

PATIENT TRANSFER SOLUTIONS



Casualty Trolleys

Introduction :

Medimek is a trusted brand in the medical fraternity over a decade. Our singular and simple philosophy is to cater to the needs of all our customers with a comprehensive range of excellent quality, user friendly hospital furniture and allied products, offered in a broad spectrum of exciting shades, aesthetic patterns and creative designs and backed by prompt after sales services.

Patient Transfer is an indispensable and round the clock activity in every medical institution, however small or large. Smooth, safe and speedy movement is vital and of paramount importance in transferring the patient in distress or trauma. Equally, effortless transfer is crucial for care givers too.

Casualty Trolley is the most suitable and dependable transfer mode for Trauma and Casualty wards and OT; especially during emergency treatment. This strong and sturdy product imparts excellent mobility; offers greater maneuverability and is equipped with varying utility accessories needed to deal with emergencies.



Mi-4076

Emergency and Recovery Trolley (Height on Hydraulic Pump and Collapsible Railings)

Specification:

- Overall approx. dimension: 2110 mm L x 730 mm W
- Stretcher top approx. dimension: 1835 mm L x 595 mm W
- Height adjustment on foot operated hydraulic pump: 660 mm to 920 mm (approx.) (without Mattress)
- Two section removable X-ray permeable stretcher top with sliding X-ray cassette holder
- Backrest raised on ratchet
- Trendelenburg and reverse trendelenburg adjustments on gas spring
- Pair of collapsible stainless steel tube railings
- Stainless steel IV pole
- Utility tray, Oxygen cylinder holder and Urine bag holder
- Buffers at four corners
- Four heavy duty castors 125 mm dia, two with brake
- Pretreated and powder coated finish
- Supplied in SKD condition

Optional Accessory:

Two section 40 mm thick PU foam mattress covered with rexine



Mi-4077

Emergency and Recovery Trolley (Height on Screw Mechanism and Collapsible Railings)

Specification:

- Overall approx. dimension: 1870 mm L x 730 mm W
 - Stretcher top approx. dimension: 1835 mm L x 595 mm W
 - Height adjustment on screw mechanism: 640 mm to 845 mm (approx.) (without Mattress)
 - Two section removable X-ray permeable stretcher top with sliding X-ray cassette holder
 - Backrest raised on ratchet
 - Trendelenburg and reverse trendelenburg adjustments on screw mechanism
 - Pair of collapsible stainless steel tube railings
 - Stainless steel IV pole
 - Oxygen cylinder holder and Urine bag holder
 - Buffers at four corners
 - Four heavy duty castors 125 mm dia, two with brake
 - Pretreated and powder coated finish
 - Supplied in SKD condition
- Optional Accessory:**
Two section 40 mm thick PU foam mattress covered with rexine



Mi-4079 A

Emergency Stretcher Trolley - SS Framework and SS Stretcher Top (Height on Screw Mechanism, SS Swing Away type Railings and Fixed Stretcher Top)

Specification:

- Overall approx. dimension: 1875 mm L x 635 mm W
- Stretcher top approx. dimension: 1735 mm L x 555 mm W
- Height adjustment on screw mechanism: 675 mm to 850 mm (approx.)
- Three section fixed stainless steel stretcher top with foot end folding
- Backrest raised on ratchet
- Trendelenburg and reverse trendelenburg adjustments on screw mechanism
- Pair of stainless steel swing away type railings
- Stainless steel IV pole, Urine bag holder
- Buffers at four corners
- Four heavy duty castors 125 mm dia, all four with brake

Mi-4079

Emergency Stretcher Trolley - MS Framework and MS Stretcher Top (Height on Screw Mechanism, MS Swing Away type Railings and Fixed Stretcher Top)

Optional accessories for Emergency Stretcher Trolley :

Mi-1037 D : Suitable 40 mm thick PU foam mattress covered with rexine for three section stretcher top

Mi-1080 : Oxygen cylinder holder (MS)

Mi-1080 A : Oxygen cylinder holder (SS)

- Accessories, colors and designs shown in the brochure illustrations may not be the part of standard product configuration.
- We reserve the right to make changes to designs, specifications and illustrations without prior notice.
- Reproduction of this brochure in any form is liable for legal action.